

Technical Publication

DICOM Conformance Statement
Digital Lightbox 2.0

Document Revision 1

February 11, 2010

2010 © Copyright BrainLAB AG



1 Conformance Statement Overview

This is a conformance statement for the BrainLAB software Digital Lightbox. One main purpose of this software is to import and visualize DICOM data.

The DICOM import part of the application is

- Browse and display of DICOM files (e.g. removable media)
- Query remote DICOM archives
- Retrieve DICOM data from archives
- Receive DICOM data from other DICOM nodes (e.g. archives or workstations) via the network
- Store DICOM data to archives

SOP Classes	User Of Service (SCU)	Provider Of Service (SCP)
Transfer		
Enhanced CT Image	Yes	Yes
Enhanced MR Image	Yes	Yes
SC Multi Frame Grayscale Byte	Yes	Yes
SC Multi Frame Grayscale Word	Yes	Yes
SC Multi Frame Single Bit	Yes	Yes
SC Multi Frame True Color	Yes	Yes
Standard CR	Yes	Yes
Standard CT	Yes	Yes
Standard Digital X-Ray Image for Image for Presentation	Yes	Yes
Standard Digital X-Ray Image for Image for Processing	Yes	Yes
Standard Grayscale Softcopy Image for Presentation State	Yes	Yes
Standard Hardcopy Color	Yes	Yes
Standard Hardcopy Grayscale	Yes	Yes
Standard Intra-oral X-Ray Image for Presentation	Yes	Yes
Standard Intra-oral X-Ray Image for Processing	Yes	Yes
Standard MG Image for Presentation	Yes	Yes
Standard MG Image for Processing	Yes	Yes
Standard MR	Yes	Yes
Standard NM	Yes	Yes
Standard NM Retired	Yes	Yes
Standard Ophthalmic 16 Bit	Yes	Yes
Standard Ophthalmic 8 Bit	Yes	Yes
Standard PET	Yes	Yes
Standard RT Image	Yes	Yes
Standard Secondary Capture	Yes	Yes
Standard US	Yes	Yes
Standard US Multi Frame	Yes	Yes
Standard US Multi Frame Retired	Yes	Yes
Standard US Retired	Yes	Yes

SOP Classes	User Of Service (SCU)	Provider Of Service (SCP)
Standard Video Endoscopic	Yes	Yes
Standard Video Microscopic	Yes	Yes
Standard Video Photographic	Yes	Yes
Standard VL Endoscopic	Yes	Yes
Standard VL Microscopic	Yes	Yes
Standard VL Photographic	Yes	Yes
Standard VL Slide Microscopic	Yes	Yes
Standard X-Ray Angio	Yes	Yes
Standard X-Ray Angio Biplane	Yes	Yes
Standard X-Ray RF	Yes	Yes
Query/Retrieve		
Patient Root Query/Retrieve Information Model - FIND	Yes	No
Patient Root Query/Retrieve Information Model - MOVE	Yes	No
Study Root Query/Retrieve Information Model - FIND	Yes	No
Study Root Query/Retrieve Information Model - MOVE	Yes	No

Table 1-1: Network services supported by Digital Lightbox

Media Storage Application Profile	Write Files (FSC or FSU)	Read Files (FSR)
Compact Disc – Recordable		
General Purpose CD–R	No	Yes

Table 1-2: Media Services supported by Digital Lightbox

2 Table Of Contents

1	Conformance Statement Overview	3
2	Table Of Contents	5
3	Introduction	7
3.1	Revision History	7
3.2	Audience	7
3.3	Remarks	7
3.4	Abbreviations	8
3.5	References	8
4	Networking	9
4.1	Implementation Model	9
4.1.1	Application Data Flow Diagram	9
4.1.2	Functional Definition of Application Entity (AE)	10
4.1.3	Sequencing Of Real World Activities	11
4.2	Application Entity Specifications	11
4.2.1	Digital Lightbox Specification	11
4.2.1.1	SOP Classes and Transfer Syntaxes	11
4.2.1.2	Association Policies	13
4.2.1.2.1	General	13
4.2.1.2.2	Number of Associations	13
4.2.1.2.3	Asynchronous Nature	14
4.2.1.2.4	Implementation Identifying Information	14
4.2.1.3	Association Initiation Policy	14
4.2.1.3.1	Activity – Find	14
4.2.1.3.2	Activity – Store	15
4.2.1.3.3	Activity – Retrieve	16
4.2.1.4	Association Acceptance Policy	16
4.2.1.4.1	Activity – Receive	17
4.3	Network Interfaces	17
4.3.1	Physical Network Interface	17
4.3.2	Additional Protocols	17
4.4	Configuration	17
4.4.1	AE Title / Presentation Address Mapping	17
4.4.1.1	Local AE Titles	18
4.4.1.2	Remote AE Title/Presentation Address Mapping	18
4.4.2	Parameters	18
5	Media Interchange	19
5.1	Implementation Model	19
5.1.1	Application Data Flow Diagram	19
5.1.2	Functional Definition of Application Entity (AE)	19
5.1.3	Sequencing Of Real World Activities	19
5.1.4	File Meta Implementation Identifying Information	19
5.2	Application Entity Specifications	20
5.2.1.1	File Meta Information for the Application Entity	20
5.2.1.2	Real-World Activities	20
5.2.1.3	Activity – Import DICOM	20
5.2.1.3.1	Media Storage Application Profiles	20
5.3	Augmented And Private Application Profiles	21

5.4	Media Configuration	21
6	Support Of Extended Character Sets	22
7	Security Profiles	23
8	Annexes	24
8.1	IOD Contents.....	24
8.1.1	Usage Of Attributes From Received IODs	24
8.1.1.1	Images	24
8.2	Data Dictionary Of Private Attributes	24
8.3	Coded Terminology And Templates.....	24
8.4	Grayscale Image Consistency	24
8.5	Standard Extended/Specialized/Private Sop Classes	24
8.6	Private Transfer Syntaxes	24
9	Indexes	25
9.1	Index Of Tables	25
9.2	Index Of Figures	25

3 Introduction

3.1 Revision History

Document Version	Date of Issue	Author	Description
3	December 7, 2009		Initial version 2.0

3.2 Audience

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

3.3 Remarks

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [1]. However, by itself it is not guaranteed to ensure the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity between BrainLAB and non-BrainLAB equipment.
- This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended. An acceptance protocol is available to validate the desired level of connectivity.
- The DICOM standard will evolve to meet the users' future requirements. BrainLAB reserves the right to make changes to its products or to discontinue its delivery.

3.4 Abbreviations

There are a variety of terms and abbreviations used in the document that are defined in the DICOM Standard. Abbreviations and terms are as follows:

AE	DICOM Application Entity
AET	Application Entity Title
CD	Compact Disk
CD-R	Compact Disk Recordable
DVD	Digital Versatile Disc
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
HD	Hard Disk
IOD	(DICOM) Information Object Definition
ISO	International Standard Organization
MOD	Magneto Optical Disk
PDU	DICOM Protocol Data Unit
Q/R	Query and Retrieve
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM server)
SOP	DICOM Service-Object Pair

3.5 References

- [1] Digital Imaging and Communications in Medicine (DICOM) 3.0, NEMA PS 3.1-3.18 – 2006

4 Networking

4.1 Implementation Model

The BrainLAB Digital Lightbox application is an implementation of:

- A Query/Retrieve SCU to query DICOM archives and to initiate a move request for the queried archive.
- A Media File Set Reader to load DICOM data from a file system.
- A Storage SCP that receives DICOM data from other DICOM archives or workstations.
- A Storage SCU that stores DICOM data to DICOM archives or workstations
- A visualization of imported DICOM data

4.1.1 Application Data Flow Diagram




Figure 4-1: The Query/Retrieve SCU application flow diagram




Figure 4-2: The Storage SCP application flow diagram




Figure 4-3: The Storage SCU application flow diagram

4.1.2 Functional Definition of Application Entity (AE)

Some communications and data transfer with remote AEs are accomplished utilizing the DICOM protocol over a network using the TCP/IP protocol stack.

- **Query and Retrieve:**
The user wants to find a certain dataset in a DICOM archive. Digital Lightbox initiates an association as a Q/R SCU negotiating all models. The find request can be performed (depending on the negotiated models) on all DICOM levels (patient, study, series or instance). A move request will be performed on series level or on instance level.
- **Storage SCP:**
With the start of the Digital Lightbox a DICOM Storage SCP is invoked. It accepts any association with a Storage SCU negotiating any of the SOP Classes listed in Table 4-3. The Storage SCP only accepts one association at a time.

Further, during a move operation a DICOM Storage SCP is invoked. Digital Lightbox accepts an association with a Storage SCU negotiating any of the SOP Classes listed in Table 4-3.

- **Storage SCU:**
The user selects the save command to send DICOM images to an archive or workstation Storage SCP. It accepts any association with a Storage SCP negotiating any of the SOP Classes listed in Table 4-3.

4.1.3 Sequencing Of Real World Activities

Digital Lightbox Query/Retrieve SCU performs a sequencing of real world activities as follows:

1. Query Archive and display received DICOM data:
 - a. Send DICOM Query/Retrieve C-FIND request.
 - b. Receive DICOM Query/Retrieve C-FIND responses.
2. User selects data to retrieve.
3. Initiate transfer of selected entities:
 - a. Send a DICOM Query/Retrieve C-MOVE request
 - b. Receive DICOM C-STORE requests with the requested SOP instances.
 - c. Receive DICOM Query/Retrieve C-MOVE responses




Figure 4-4: Sequencing of Query/Retrieve SCU and Storage SCP

4.2 Application Entity Specifications

4.2.1 Digital Lightbox Specification

4.2.1.1 SOP Classes and Transfer Syntaxes

Digital Lightbox sends or receives a C-ECHO request in order to test the connection to a remote AE. It provides standard conformance to the following DICOM V3.0 SOP Classes:

SOP Class Name	SOP Class UID	SCU	SCP
Verification SOP Class	1.2.840.10008.1.1	Yes	Yes

Table 4-1: Supported Verification SOP Classes

Digital Lightbox is able to query a remote archive. It provides Standard Conformance to the following DICOM V3.0 SOP Classes:

SOP Class Name	SOP Class UID	SCU	SCP
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1	Yes	No
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2	Yes	No
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1	Yes	No
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2	Yes	No

Table 4-2: Supported Query/Retrieve SOP Classes

Digital Lightbox imports and exports DICOM image data. It provides Standard Conformance to the following DICOM V3.0 SOP Classes:

SOP Class Name	SOP Class UID	SCU	SCP
Enhanced CT Image	1.2.840.10008.5.1.4.1.1.2.1	Yes	Yes
Enhanced MR Image	1.2.840.10008.5.1.4.1.1.4.1	Yes	Yes
SC Multi Frame Grayscale Byte	1.2.840.10008.5.1.4.1.1.7.2	Yes	Yes
SC Multi Frame Grayscale Word	1.2.840.10008.5.1.4.1.1.7.3	Yes	Yes
SC Multi Frame Single Bit	1.2.840.10008.5.1.4.1.1.7.1	Yes	Yes
SC Multi Frame True Color	1.2.840.10008.5.1.4.1.1.7.4	Yes	Yes
Standard CR	1.2.840.10008.5.1.4.1.1.1	Yes	Yes
Standard CT	1.2.840.10008.5.1.4.1.1.2	Yes	Yes
Standard Digital X-Ray Image for Image for Presentation	1.2.840.10008.5.1.4.1.1.1.1	Yes	Yes
Standard Digital X-Ray Image for Image for Processing	1.2.840.10008.5.1.4.1.1.1.1.1	Yes	Yes
Standard Grayscale Softcopy Image for Presentation State	1.2.840.10008.5.1.4.1.1.11.1	Yes	Yes
Standard Hardcopy Color	1.2.840.10008.5.1.1.30	Yes	Yes
Standard Hardcopy Grayscale	1.2.840.10008.5.1.1.29	Yes	Yes
Standard Intra-oral X-Ray Image for Presentation	1.2.840.10008.5.1.4.1.1.1.3	Yes	Yes
Standard Intra-oral X-Ray Image for Processing	1.2.840.10008.5.1.4.1.1.1.3.1	Yes	Yes
Standard MG Image for Presentation	1.2.840.10008.5.1.4.1.1.1.2	Yes	Yes
Standard MG Image for Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Yes	Yes
Standard MR	1.2.840.10008.5.1.4.1.1.4	Yes	Yes
Standard NM	1.2.840.10008.5.1.4.1.1.20	Yes	Yes
Standard NM Retired	1.2.840.10008.5.1.4.1.1.5	Yes	Yes
Standard Ophthalmic 16 Bit	1.2.840.10008.5.1.4.1.1.77.1.5.2	Yes	Yes
Standard Ophthalmic 8 Bit	1.2.840.10008.5.1.4.1.1.77.1.5.1	Yes	Yes
Standard PET	1.2.840.10008.5.1.4.1.1.128	Yes	Yes
Standard RT Image	1.2.840.10008.5.1.4.1.1.481.1	Yes	Yes
Standard Secondary Capture	1.2.840.10008.5.1.4.1.1.7	Yes	Yes
Standard US	1.2.840.10008.5.1.4.1.1.6.1	Yes	Yes
Standard US Multi Frame	1.2.840.10008.5.1.4.1.1.3.1	Yes	Yes

SOP Class Name	SOP Class UID	SCU	SCP
Standard US Multi Frame Retired	1.2.840.10008.5.1.4.1.1.3	Yes	Yes
Standard US Retired	1.2.840.10008.5.1.4.1.1.6	Yes	Yes
Standard Video Endoscopic	1.2.840.10008.5.1.4.1.1.77.1.1.1	Yes	Yes
Standard Video Microscopic	1.2.840.10008.5.1.4.1.1.77.1.2.1	Yes	Yes
Standard Video Photographic	1.2.840.10008.5.1.4.1.1.77.1.4.1	Yes	Yes
Standard VL Endoscopic	1.2.840.10008.5.1.4.1.1.77.1.1	Yes	Yes
Standard VL Microscopic	1.2.840.10008.5.1.4.1.1.77.1.2	Yes	Yes
Standard VL Photographic	1.2.840.10008.5.1.4.1.1.77.1.4	Yes	Yes
Standard VL Slide Microscopic	1.2.840.10008.5.1.4.1.1.77.1.3	Yes	Yes
Standard X-Ray Angio	1.2.840.10008.5.1.4.1.1.12.1	Yes	Yes
Standard X-Ray Angio Biplane	1.2.840.10008.5.1.4.1.1.12.3	Yes	Yes
Standard X-Ray RF	1.2.840.10008.5.1.4.1.1.12.2	Yes	Yes

Table 4-3: Supported Storage SOP Classes

Digital Lightbox supports the following transfer syntaxes. In an association negotiation the syntaxes are proposed in the order of appearance in the following tables.

Transfer Syntax Name	Transfer Syntax UID	SCU	SCP
DICOM Implicit VR Little Endian	1.2.840.10008.1.2	Yes	Yes
DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1	Yes	Yes
DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2	Yes	Yes

Table 4-4: Supported Standard Transfer Syntaxes

Transfer Syntax Name	Transfer Syntax UID	SCU	SCP
JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1]): Default Transfer Syntax for Lossless JPEG Image Compression	1.2.840.10008.1.2.4.70	Yes	Yes

Table 4-5: Supported Compression Transfer Syntaxes

4.2.1.2 Association Policies

4.2.1.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.1.2.2 Number of Associations

Maximum number of simultaneous Associations (Initiator)	2 (configurable)
Maximum number of simultaneous Associations (Acceptor)	unrestricted

4.2.1.2.3 Asynchronous Nature

The Digital Lightbox does not support asynchronous communication (multiple outstanding transactions over a single association).

Maximum number of outstanding asynchronous transactions	1
---	---

4.2.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Implementation Class UID	1.2.276.0.20.1.1.12.2.0.0
Implementation Version Name	DIGITALLIGHTBOX2

4.2.1.3 Association Initiation Policy

Digital Lightbox initiates an association in these cases:

1. Find: The user tries to find a specific entity in a remote DICOM archive.
2. Store: The user wants to send a specific entity to a remote DICOM archive.
3. Retrieve: The user wants to retrieve a specific entity from the remote DICOM archive.

4.2.1.3.1 Activity – Find

4.2.1.3.1.1 Description and Sequencing of Activities

A DICOM Query/Retrieve C-FIND request is performed when the user queries the remote DICOM archive for patients, studies, series or instances.

4.2.1.3.1.2 Proposed Presentation Contexts

Presentation Context Table			
Abstract Syntax	Transfer Syntax	Role	Ext. Neg.
All SCU SOP Classes listed in Table 4-1 and all FIND SCU SOP Classes as listed in Table 4-2	All SCU Transfer Syntaxes as listed in Table 4-4	SCU	None
		SCU	None
		SCU	None

Table 4-6: Proposed Presentation Contexts for Activity Find.

4.2.1.3.1.3 SOP Specific Conformance

Digital Lightbox provides standard conformance to the DICOM Verification Service Class and to the DICOM Query/Retrieve FIND SOP Classes. No extended negotiation is implemented.

A C-FIND can be executed on all levels. On a lower level, all key attributes of the higher levels are included in the query (i.e. no hierarchical queries are performed). C-FIND's on IMAGE level can be disabled via a configuration option.

A C-CANCEL will be sent if the user aborts the search process or the Digital Lightbox is shut down. If the remote server does not react to a C-CANCEL, the association is aborted.

Attribute Name	Tag	Types of Matching
Patient Level (Patient Root Q/R Information Model) Study Level (Study Root Q/R Information Model)		
Patient's Name	(0010,0010)	S, *, U
Patient ID	(0010,0020)	S, *, U
Patient's Birth Date	(0010,0030)	U
Patient's Sex	(0010,0040)	U
Study Level		
Study Date	(0008,0020)	S, R, U
Accession Number	(0008,0050)	S, *, U
Study Instance UID	(0020,000D)	UNIQUE
Study Time	(0008,0030)	U
Study Description	(0008,1030)	U
Series Level		
Series Instance UID	(0020,000E)	UNIQUE
Series Number	(0020,0011)	U
Series Description	(0008,103E)	U
Modality	(0008,0060)	U
No. of Series Rel. Instances	(0020,1209)	U
Instance Level		
SOP Instance UID	(0008,0018)	UNIQUE
Instance Number	(0020,0013)	U

Table 4-7: Patient Root and Study Root Request Identifier for FIND-SCU

The types of Matching supported by the C-FIND SCU:

- An "S" indicates the identifier attribute uses Single Value Matching.
- An "R" indicates Range Matching.
- A "*" indicates wildcard matching.
- A 'U' indicates Universal Matching (i.e. the attributes value is empty and shall only be returned, not used as a filter key).
- "UNIQUE" indicates that this is the Unique Key for that query level, in which case Universal Matching or Single Value Matching is used depending on the query level.

4.2.1.3.2 Activity – Store

4.2.1.3.2.1 Associated Real-World Activity

To send DICOM data a storage request is performed to a remote Storage SCP. The remote Storage SCP must be one of the AETs known by Digital Lightbox.

4.2.1.3.2.2 Proposed Presentation Contexts

Presentation Context Table			
Abstract Syntax	Transfer Syntax	Role	Ext. Neg
All SCU SOP Classes as listed in Table 4-1 and Table 4-3	All SCU Transfer Syntaxes as listed in Table 4-4 and Table 4-5	SCU	None
		SCU	None

Table 4-8: Proposed Presentation Contexts for Activity Store.

4.2.1.3.2.3 SOP Specific Conformance

Digital Lightbox provides standard conformance to the DICOM Verification Service Class and to the DICOM Storage SOP Classes. No extended negotiation is implemented.

4.2.1.3.3 Activity – Retrieve

4.2.1.3.3.1 Description and Sequencing of Activities

On user selection of a specific DICOM study or series (depends on the available SOP Class for Query/Retrieve), a move request is performed. The storage target for receiving the DICOM data (the AET with which the move-request is equipped) is the Digital Lightbox application itself.

The Move operation is only performed on SERIES or IMAGE level. Image level retrieval is used for preview image loading and can be disabled via a configuration setting.

A C-CANCEL will be sent if the user selects a different image series or if the Digital Lightbox is shut down. If the remote server does not react to a C-CANCEL, the association is aborted.

The Move operation only can be invoked after a Find operation. See chapter 4.1.3 for a detailed sequencing diagram.

4.2.1.3.3.2 Proposed Presentation Contexts

Presentation Context Table			
Abstract Syntax	Transfer Syntax	Role	Ext. Neg.
All SCU SOP Classes listed in Table 4-1 and all SCU MOVE SOP Classes as listed in Table 4-2	All SCU Transfer Syntaxes as listed in Table 4-4	SCU	None
		SCU	None
		SCU	None

Table 4-9: Proposed Presentation Contexts for Activity Move.

4.2.1.3.3.3 SOP Specific Conformance

Digital Lightbox provides standard conformance to the DICOM Verification Service Class and to the DICOM Query/Retrieve MOVE SOP Classes. No extended negotiation is implemented.

4.2.1.4 Association Acceptance Policy

Digital Lightbox accepts an association in this case:

1. Receive:
 - The user wants to retrieve a specific entity from the remote DICOM archive.
 - The Digital Lightbox responds to storage requests.

4.2.1.4.1 Activity – Receive

4.2.1.4.1.1 Associated Real-World Activity

As DICOM storage instances are received they are saved to the local file system. If the received instance is a duplicate of a previously received instance, the new instance will not be stored.

4.2.1.4.1.2 Proposed Presentation Contexts

Presentation Context Table			
Abstract Syntax	Transfer Syntax	Role	Ext. Neg
All SCP SOP Classes as listed in Table 4-1 and Table 4-3	All SCP Transfer Syntaxes as listed in Table 4-4 and Table 4-5	SCP	None
		SCP	None
		SCP	None

Table 4-10: Storage SCP Presentation Contexts.

4.2.1.4.1.3 SOP Specific Conformance

The Digital Lightbox provides standard conformance to the DICOM Verification Service Class and to the DICOM Storage SOP Classes. No extended negotiation is implemented.

4.2.1.4.1.4 Presentation Context Acceptance Criterion

The Digital Lightbox accepts multiple presentation contexts containing the same abstract syntax.

4.2.1.4.1.5 Transfer Syntax Selection Policy

The first Transfer Syntax encountered in the configuration file, which matches a Transfer Syntax offered for a given Presentation Context, will be selected as the accepted Transfer Syntax for that Presentation Context.

4.3 Network Interfaces

4.3.1 Physical Network Interface

Digital Lightbox supports the DICOM upper layer using TCP/IP. Digital Lightbox is indifferent to the physical medium over which TCP/IP executes. It inherits this from the operating system upon which it executes.

4.3.2 Additional Protocols

The usage of DNS and DHCP is possible and is based on the network configuration of the operating system upon which Digital Lightbox executes.

4.4 Configuration

4.4.1 AE Title / Presentation Address Mapping

Configuration of remote and local DICOM nodes can be performed with the graphical user interface of Digital Lightbox settings configuration program.

Digital Lightbox can configure several nodes representing remote Q/R Servers. On the corresponding settings page, application-wide global parameter and node-specific parameters can be entered.

4.4.1.1 Local AE Titles

The Local AET and the Move Destination AET for the Query/Retrieve SCU can be configured for each remote Q/R Server.

The AET of the Storage SCP does not need to be configured, since the Digital Lightbox accepts any called AET. It is recommended to use LIGHTBOX_STORE.

Application Entity	Default AE Title	Default TCP/IP Port
Digital Lightbox (Query/Retrieve)	LIGHTBOX_QR	104
Digital Lightbox (Storage SCP only)	LIGHTBOX_STORE	104

4.4.1.2 Remote AE Title/Presentation Address Mapping

In Digital Lightbox you can specify several nodes for import. The IP address/hostname, AET and listening port may be configured for each DICOM network node separately within the Digital Lightbox settings configuration program's graphical user interface.

4.4.2 Parameters

Additional a timeout may be specified for each DICOM network archive separately.

Parameter	Configurable	Default Value
Timeout	Yes	30
Maximum PDU Size	No	64234

5 Media Interchange

Digital Lightbox supports DICOM media interchange for import of DICOM data:

- For import Digital Lightbox supports media interchange application profiles. To reflect this, the support for the Standard General Purpose CD-R Interchange is added to provide the supported SOP Classes. Nevertheless Digital Lightbox is able to import DICOM files even without the existence of any DICOMDIR by scanning a given file system located on any media (e.g. HD, MOD, CD, DVD, Tapes, USB Drive) for any kind of DICOM files.

5.1 Implementation Model

5.1.1 Application Data Flow Diagram

With Digital Lightbox the user may browse DICOM file sets and import selected entities.




Figure 5-1: The media interchange application flow diagram

5.1.2 Functional Definition of Application Entity (AE)

Some communications and data transfer with remote AEs are accomplished utilizing the file system provided by the operating system upon which Digital Lightbox executes.

- File Set Reader:
Digital Lightbox loads DICOM data from the file. The reader supports the same SOP classes as the Storage SCP (see Table 4-3).

5.1.3 Sequencing Of Real World Activities

Not necessary.

5.1.4 File Meta Implementation Identifying Information

Digital Lightbox provides the same information as in chapter 4.2.1.2.4.

5.2 Application Entity Specifications

Digital Lightbox supports the following Media Interchange Profiles:

AE Related Application Profiles, Real-World Activities, and Roles			
Supported APs	Real World Activity	Roles	SC Option
STD-GEN-CD	Import DICOM	FSR	Interchange

Table 5-1: Supported Media Interchange Profiles.

5.2.1.1 File Meta Information for the Application Entity

The Source Application Entity Title included in the File Meta Header is configurable. It is the same as the local AET of the network configuration (see section 4.4.1.1)

5.2.1.2 Real-World Activities

5.2.1.3 Activity – Import DICOM

Digital Lightbox acts as an FSR using the Interchange option

- When requested to provide a directory listing it will read the File-set and display the DICOM-DIR directory entries for all SOP Instances in the File-set.
- When requested to import the selected entries from directory listing, only those SOP Instances are loaded that correspond to the Application Profile STD-GEN-CD.
- For the list of Application Profiles invoking this AE see Table 5-1. The supported SOP Classes see Table 4-3.

5.2.1.3.1 Media Storage Application Profiles

Digital Lightbox supports the STD-GEN-CD Application Profile.

5.2.1.3.1.1 Options

Supported transfer syntaxes for the media profiles:

Transfer Syntax Name	Transfer Syntax UID
DICOM Implicit VR Little Endian	1.2.840.10008.1.2
DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1
DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2
JPEG Baseline (Process 1): Default Transfer Syntax for Lossy JPEG 8 Bit Image Compression	1.2.840.10008.1.2.4.50
JPEG Extended (Process 2 & 4): Default Transfer Syntax for Lossy JPEG 12 Bit Image Compression (Process 4 only)	1.2.840.10008.1.2.4.51
JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1]): Default Transfer Syntax for Lossless JPEG Image Compression	1.2.840.10008.1.2.4.70
RLE Lossless	1.2.840.10008.1.2.5

Table 5-2: Supported Media Profile Transfer Syntaxes

The Offline-Media Application Entity supports the SOP Classes and Transfer Syntaxes listed in the Table below:

SOP Classes	Transfer Syntaxes
All SCP SOP Classes listed in Table 4-3	All SCP Transfer Syntaxes listed in Table 5-2

For further information see section 8.1.1 on acceptance of SOP Instances (i.e. whether Digital Lightbox is able to import and convert the DICOM data).

5.3 Augmented And Private Application Profiles

Digital Lightbox does not support any augmented or private application profiles.

5.4 Media Configuration

Digital Lightbox uses the local AET configured for the Store SCP network service as source AET for the DICOM files.

6 Support Of Extended Character Sets

Digital Lightbox supports the

- ISO_IR 100 (ISO 8859-1:1987 Latin Alphabet No. 1 supplementary set)

7 Security Profiles

No security profiles are supported.

8 Annexes

8.1 IOD Contents

8.1.1 Usage Of Attributes From Received IODs

This section describes the requirements on the DICOM data, which can be displayed.

8.1.1.1 Images

Digital Lightbox accepts all images of the SOP Classes in Table 4-3. Though, there are some restrictions and special conversions:

- Images with an attribute (0028,0030) Pixel Spacing containing different values for x and y distance¹ will be ignored.

8.2 Data Dictionary Of Private Attributes

None supported.

8.3 Coded Terminology And Templates

None supported.

8.4 Grayscale Image Consistency

Not supported.

8.5 Standard Extended/Specialized/Private Sop Classes

None supported.

8.6 Private Transfer Syntaxes

None supported.

¹ To be more precise: If the difference between both values is greater than 0.001 mm.

9 Indexes

9.1 Index Of Tables

Table 1-1: Network services supported by Digital Lightbox	4
Table 1-2: Media Services supported by Digital Lightbox	4
Table 4-1: Supported Verification SOP Classes	11
Table 4-2: Supported Query/Retrieve SOP Classes.....	12
Table 4-3: Supported Storage SOP Classes	13
Table 4-4: Supported Transfer Syntaxes.....	13
Table 4-5: Proposed Presentation Contexts for Activity Find.....	14
Table 4-6: Patient Root and Study Root Request Identifier for FIND-SCU	15
Table 4-7: Proposed Presentation Contexts for Activity Store.....	16
Table 4-8: Proposed Presentation Contexts for Activity Move.....	16
Table 4-9: Storage SCP Presentation Contexts.....	17
Table 5-1: Supported Media Interchange Profiles.....	20
Table 5-2: Supported Media Profile Transfer Syntaxes	20

9.2 Index Of Figures

Figure 4-1: The Query/Retrieve SCU application flow diagram	9
Figure 4-2: The Storage SCP application flow diagram	10
Figure 4-3: The Storage SCU application flow diagram.....	10
Figure 4-4: Sequencing of Query/Retrieve SCU and Storage SCP	11
Figure 5-1: The media interchange application flow diagram.....	19