

# WHAT'S NEW IN VERSION 2.4

---



## TOUCH SCREEN MODE

Touch screen mode allows you to plan on Windows 7 & 8 touch displays.

Using multi-touch gestures, zoom and pan the image, select, rotate and adjust any object, change implant size, and perform measurements.



## QENTRY INTEGRATION

Template from any internet-connected Mac or PC, using TraumaCad's integration with Qentry image sharing. Qentry provides a secure online environment to view, access, and share patient images as well as jpeg, PDF and doc files.



## NEW KNEE PACKAGE

New wizards for faster pre-operative AP resection line measurement and post-operative results evaluation.



## ADVANCED SPINE PACKAGE

New and enhanced measurements for quick and accurate spine pre-operative planning on AP and LAT images.

