



WE HAVE
SOME
EXCITING
NEWS!



A NEW
NAME!

FROM VOYANT
HEALTH TO
BRAINLAB LTD.

A NEW
WEBSITE!

CHECK OUT
TRAUMACAD.COM

NEW
RELEASES

TRAUMACAD 2.5
TRAUMACAD
MOBILE 2.0

NEW FEATURES

WE KNOW YOUR TIME IS VALUABLE, SO WE HAVE WORKED HARD TO MAKE THE PLANNING PROCESS EASIER BY ADDING NEW AUTOMATIC FEATURES AND ENHANCEMENTS THAT YOU HAVE ASKED FOR:

- **AUTO KNEE** automatically places resection lines, positions templates, estimates size, and simulates the alignment.
- **AUTO HIP MEASUREMENTS** Cup Version tool is now automatically positioned.
- **SAVED IMAGE REPORT** define what information will be displayed on the planned image, and customize the font and background.
- **TEMPLATE PROPERTIES WINDOW** has now moved to the left panel. At any point, you can bring the floating window back from the Settings menu.
- **NEW PLANNING PROCEDURES** for TraumaCad Mobile: Knee, Upper Limb and Foot & Ankle, including a double image planning functionality.